




## P-Mill Grinding for Laboratory

### Horizontal Mill for Laboratory



P-Mill is used in research & development for dispersion & wet grinding

- ❖ Compact size
- ❖ Easy to operate
- ❖ High efficiency
- ❖ Manual/Pump feed

#### Option

- Explosion proof type of motor and control panel
- Temperature controller
- Gear pump set

Model	Shell Volume (liters)	Vessel Capacity (liters)	Out-let l/h	Mill Motor (HP)	A Width (mm)	B Length (mm)	C Height (mm)	Weight (Kg)
P-Mill 0.5 / 0.5 Ex	1.1	0.5	3 - 5	2	500	750	630	90 / 101
P-Mill 1 / 1 Ex	1.1	1.0	5 - 10	2	500	750	630	107 / 110
P-Mill 2 / 2 Ex	2.2	2.0	5 - 20	3	500	850	630	110 / 119

P.General Products : 234 Moo.5,Prakasa Road.,Prakasa,Muang,Samutprakarn 10280 Thailand.

Tel: (662) 728-5246-8 Fax: (662) 728-5249

E-mail: info@p-general.com